


MEASURING TEMPERATURE

4. Why do thermometers need to be handled carefully?
- made of glass, break easily
 - may have mercury inside them if they are old
 - are round and can roll (off the table)
5. a) Suppose you are working at location “A”. Draw (on the diagram) how a thermometer should be placed on a lab bench when you are not using it.
- b) Why would working with a thermometer at “B” be a bad idea?



- no way to store it perpendicular to the table...will always be able to roll
- a potentially high-traffic zone because close to the sink, may have many people brushing by you as you work

MEASURING MASS

6. Label the following with the names and functions of the equipment. (If you are unsure, ask!)

triple beam balance	digital (electronic) balance	weigh boat	scoopula
			
weighing objects	weighing objects	hold substance being weighed	transfer small quantities of solids (especially when weighing)
_____	_____	_____	_____

Using a Digital Balance

There are a number of steps to using a digital/electronic balance correctly. Practice weighing out 5.00 grams of the substance provided. Then, if needed, add your own notes and/or diagrams to help you remember what the steps mean.

- 1) Turn on the balance.
- 2) Put a weigh boat onto the balance. Press Tare.
- 3) Use a scoopula to weigh out your substance directly onto the weigh boat.
 - a. Take care not to drop any outside the weigh boat.
 - b. Slow down with your additions when you get closer and closer to your desired amount. If you go over, **do not return the substance to the original bottle**: follow your teacher’s instructions to discard the excess.
- 4) To weigh out a different substance, gently brush off anything that is on the surface of the balance, and repeat steps 2 and 3.

Practice Questions

7. What does it mean to tare a balance? Why is this important when weighing things?

Taring means to set the 'current weight' (of the atmospheric gases on top of the scale) as equal to 0. This is important so that we can avoid including the weight of the container in the readout.

8. You are being asked to measure out 10 grams of sugar for your experiment. You follow all the steps and your scale reads 10.00 grams. But after you take off your weigh boat to pour the sugar into the beaker, you see that there are some sugar crystals left on the scale. Will this affect your experiment? If so, how, and what should you do?

Yes this will affect the experiment. We will end up with less than 10.00g of sugar.

Consult the teacher. There may be a way to figure out how much was on the scale so that we can make up the difference and add that to the sugar we have already weighed. If not, we will have to dispose of the sugar and restart, remembering to tare this time.

9. Why do we never return excess chemicals to the original bottles?

Risk of contamination.

Returning excess chemical could also cause unexpected reactions to take place, if the chemical we are returning was labelled incorrectly or was contaminated.

10. The digital balance is not usually used to weigh liquids. Why do you think this is?






Electronics and water do not mix.

MEASURING VOLUMES

Volume is the amount of space something takes up. You are probably familiar with measuring volumes of things at home in your kitchen (e.g. 1 cup of sugar; 150 mL of milk).

- Volume of a solid is measured in cm³ or centimeters cubed.
- Volume of a liquid is measured in mL or milliliters.
- Note: 1 cm³=1mL. (Especially important when using the water displacement method to measure the volume of solids.)

11. Measuring the volume of a liquid can be done using laboratory equipment with markings on the sides. Examples are below. Label these pieces of equipment (not drawn to scale).


dropper	graduated cylinder	pipette	buret	beaker
				

12. Match the following with their functions.

- a. Pipette
 - b. Graduated cylinder
 - c. Beaker
 - d. Buret
 - e. Dropper
- Measuring large volumes (e.g. 40 mL, 120 mL)
Measuring small volumes (e.g. 5 mL, 11 mL)
Measuring volumes very precisely (e.g. 10.3 mL, 34.5 mL)
Transferring very small volumes (not very precise)


In a science lab, we want to make sure that we are measuring *exactly* the amount that is asked of us. If we don't, experiments may not run the way we expect, or we may waste a lot of (sometimes extremely very expensive!) chemicals.

Sometimes, this can be difficult in graduated cylinders, burets, and pipettes. In narrow tubes, liquids will form a **meniscus**: the curved upper surface of the liquid. A meniscus forms because of attraction between molecules (when this happens between water molecules, this is referred to as water tension.)



To accurately measure volumes, you must:

- 1) Make sure your eye is level with the liquid level.
- 2) Read the volume that lines up with the bottom of the meniscus.
- 3) Record your volume accurately: there should be one estimated digit.



A graduated cylinder with 36.5 mL of water. Make sure your eyes are lined up with the level of the liquid when recording.

Practice Questions

13. Record the volumes in these graduated cylinders. Units are in mL.

1) 24.8mL

2) 9.0mL

3) 40.0mL

4) 20.0mL

14. Using one of the graduated cylinders available, practice precisely measuring out 20.0mL of tap water.

Water Displacement Method for Solids

The volume of solids can be calculated (using math formulas) or measured. Here, we will focus on measurement using the water displacement method, which is most useful for irregularly shaped solids but can be used for liquids too.

- 1) Select a container (e.g. beaker, graduated cylinder) the object will fit into.
- 2) Pour water into the empty container until it is about half full.
- 3) Record the volume of water in the container.
- 4) Carefully place the object into the container, and make sure it is fully submerged.
- 5) Record the new volume of the water plus the object.
- 6) Calculate the volume of the solid using the formula:

$$\text{volume}_{\text{solid}} = \text{volume}_{\text{final}} - \text{volume}_{\text{initial}}$$

- 7) Record the volume with proper units. Remember: $1\text{mL}=1\text{cm}^3$. mL is used for liquids, and cm³ for solids.

15. Practice measuring the volumes of the provided solids, using the water displacement method. Remember to clean up each station after you are done (wipe up spilled water with paper towel).

Station	Description of Object	Volume (with units!)
1	marble	2.6 cm ³
2	ping pong ball	32 cm ³
3	eraser	13.7 cm ³
4	stress ball	174 cm ³

16. Why is it important that the object be fully submerged before you record the total volume (step 4)?
or else your calculated volume will not be accurate (it will be less than the actual volume)

17. Draw diagrams to illustrate each of the steps in determining the volume of a solid. Use the back of the page if necessary.